
# REF M46 MICROGEN® CAMPYLOBACTER

GB

## INTENDED USE

Microgen® Campylobacter is a rapid latex agglutination test intended for confirmatory identification of enteropathogenic, thermophilic campylobacters cultured on selective solid media from faecal samples from patients with suspected bacterial enteritis. The kit is intended for professional laboratory use only.

## PRINCIPLE OF THE TEST

Latex particles are coated with rabbit immunoglobulins raised against antigen preparations from selected *Campylobacter jejuni* serotypes. When the sensitised latex particles are mixed with a solution containing enteropathogenic Campylobacter antigens, a sensitive and specific immunochemical reaction takes place causing the finely dispersed latex particles to agglutinate into aggregates that are easily visible to the unaided eye.

CONT

## KIT PRESENTATION

REAG	TEST
------	------

M46a Test Latex Reagent: 2.5mL

Latex particles coated with rabbit antibodies to Campylobacter antigens. Preserved with 0.099% sodium azide. (Blue cap)

REAG	CONTROL
------	---------

M46b Control Latex Reagent: 2.5mL

Latex particles coated with non-specific rabbit immunoglobulins. Preserved with 0.099% sodium azide. (Yellow cap)

CONTROL	+
---------	---

M46c Positive Control: 1.0mL

Suspension of inactivated Campylobacter antigens reactive with Test Latex Reagent and non-reactive with Control Latex Reagent. Preserved with 0.099% sodium azide. (Black cap)

SAMP	DIL
------	-----

M46d Sample Diluent: 50mL

0.85% Isotonic Saline preserved with 0.099% sodium azide.(white cap)

Instructions for Use  
Disposable agglutination slides  
Disposable mixing sticks

## Additional Requirements:

- Bacteriological loops
- Campylobacter semi-selective solid agar culture plates (see below)
- Disposable plastic specimen tubes (capacity 1mL)
- Gas jars for maintaining microaerophilic culture conditions
- 42°C incubator (if not available, a 37°C incubator is acceptable)

## WARNINGS AND PRECAUTIONS

## Safety:

1. The reagents supplied in this kit are for *in vitro* diagnostic use only
2. Sodium azide, which is used as a preservative in the kit reagents can react with lead or copper plumbing to form potentially explosive metal azides. Dispose by flushing with a large volume of water to prevent azide build-up.
3. Appropriate precautions should be taken when handling or disposing of potential pathogens. Decontamination of infectious material can be achieved with sodium hypochlorite at a final concentration of 3% for 30 minutes. Liquid waste containing acid must be neutralised before treatment.
4. The positive control has been inactivated during the manufacturing process. However, it should be handled as though potentially infectious.

## Procedural:

1. Microgen® Campylobacter should be used according to the kit instructions.
2. Allow all reagents to reach room temperature before use.
3. Do not dilute any of the kit reagents
4. Do not intermix reagents from different batches of kits.
5. Do not freeze any of the kit reagents
6. Do not allow the latex reagent dropper to touch the positive control or bacterial samples.
7. Ensure the agglutination slide is clean and dry prior to use.

## STORAGE AND SHELF LIFE

Microgen® Campylobacter should be stored at 2-8°C when not in use. The kit should not be used after the expiry date printed on the carton label.

## SPECIMENS

Faecal samples should be inoculated on to blood-free selective agar plates (Modified CCDA Preston, see Appendix) at a concentration of 0.2-0.3g sample per plate. Plates should be incubated in a microaerophilic atmosphere at 42°C for 48 hours. Colonies with morphology resembling Campylobacter are removed for testing with Microgen® Campylobacter (see below)

## PROCEDURE

### Quality Control:

The following check with the Positive Control should be performed each time the kit is used to confirm that the reagents are functioning correctly:

1. A single 50µL drop of Positive Control should be dispensed on to two adjacent areas on the test slide.
2. These should be tested with the Test and Control Latex Reagents as described in "Test Procedure" below.
3. Deterioration of a reagent should be suspected if:
  - There is no reaction between the Test Latex Reagent and the Positive Control OR the reaction shows a significant loss of strength with time.
  - The Control Latex Reagent reacts with the Positive Control
  - A latex reagent becomes discoloured or forms lumps which do not disperse with gentle shaking.

### Test Procedure:

1. Bring all reagents to room temperature. Gently shake the latex reagents to ensure a homogeneous suspension.
2. Dispense 50µL of Sample Diluent on to each of two ovals of the agglutination slide.
3. Use an inoculating loop to remove several colonies with Campylobacter-like morphology. If evidence of microbial growth is sparse, take a broad sweep of the agar surface. Mix the bacteria into each of the two drops of Sample Diluent on the slide to form even suspensions.
4. Add 1 drop (50µL) of Control Latex reagent to one of the bacterial suspensions on the slide. Similarly dispense 1 drop (50µL) of Test Latex reagent to the other bacterial suspension.
5. Mix the bacterial suspensions with the latex reagents using a mixing stick, **starting with the Control Latex reagent**. Spread the mixtures to the edges of the oval areas.
6. Gently rock the slide, to keep the fluid suspensions in constant movement, for 2 minutes. Observe for agglutination.
7. Read the test results (see INTERPRETATION below)
8. Discard the used mixing sticks and slides into a suitable disinfectant.

### INTERPRETATION

An agglutination reaction is indicated by visible aggregation of the latex particles. The strength of the reaction may vary, and can be assessed according to the following guidelines.

- + reaction - fine, but readily discernible granularity against a milky background.
- ++ reaction - coarse granularity against a milky background.
- +++ reaction - heavy clumping of particles around the periphery of the test oval, against a clear background.

Microgen® Campylobacter results should be interpreted as follows:

Reaction with Test Latex	Reaction with Control Latex	Interpretation
+	-	Campylobacters present
-	-	Campylobacters not present in sufficient numbers to be detected by the test
+	+	Non-specific agglutination. Inconclusive result*

\*A specimen that causes the Control Latex reagent to agglutinate cannot be tested by Microgen® Campylobacter.

### LIMITATIONS OF USE

1. Results should be interpreted by the clinician in the context of all available clinical and laboratory information.
2. Non-specific agglutination of the Control Latex reagent does not preclude the presence of Campylobacters but the result has to be reported as inconclusive. The sample should be tested by an alternative method.
3. Very low numbers of Campylobacters may result in a negative test result. 48-hour cultures should be used to maximise growth of the bacteria.

### PERFORMANCE CHARACTERISTICS

Microgen® Campylobacter has been evaluated as a culture confirmation test at two independent external centres and in-house. In total, 527 faecal specimens were cultured on selective agar plates for 48 hours and colonies tested by Microgen® Campylobacter.

		Microgen® Campylobacter		Total
		Positive	Negative	
Reference Methods	Positive	143	2	145
	Negative	1	381	382
Total		144	383	527

Sensitivity: 143/145 = 98.6%  
 Specificity: 381/382 = 99.7%

Efficiency: 524/527 = 99.4%

The specificity of Microgen® Campylobacter has been confirmed by testing a wide variety of cultured microorganisms:

Organism	Number of isolates	Microgen® Campylobacter		
		Pos	Neg	Non-specific reaction
<b>Thermophilic Campylobacters:</b>				
<i>C. jejuni</i>	37	37	0	0
<i>C. jejuni</i> subsp. <i>doylei</i>	3	3	0	0
<i>C. coli</i>	4	4	0	0
<i>C. upsaliensis</i>	13	13	0	0
<i>C. laridis</i>	1	1	0	0
<i>C. fetus</i>	2	2	0	0
<b>Non-thermophilic campylobacters:</b>				
<i>C. concisus</i>	5	0	3	2 (weak)
<i>C. hyointestinalis</i>	3	0	3	0
<i>C. sputorum faecalis</i>	1	0	0	1 (weak)
<i>C. curvis</i>	2	0	2	0
<b>Closely related bacteria:</b>				
<i>Helicobacter pylori</i>	12	2 (weak)	10	0
<i>Helicobacter fennelliae</i>	5	0	5	0
<i>Helicobacter cinaedi</i>	4	2	2	0
<i>Arcobacter butzleri</i>	4	0	4	0
<i>Arcobacter skirrowii</i>	1	0	1	0
<b>Other non-related bacteria:</b>				
<i>Providencia stuartii</i>	1	0	1	0
<i>Providencia acalifaciens</i>	1	0	1	0
<i>Escherichia coli</i>	5	0	5	0
<i>Pseudomonas aeruginosa</i>	2	0	2	0
<i>Citrobacter spp</i>	9	0	9	0
<i>Citrobacter freundii</i>	1	0	1	0
♣ Streptococcus	2	0	2	0
β Streptococcus	1	0	1	0
<i>Klebsiella pneumoniae</i>	1	0	1	0
<i>Klebsiella edwardsii</i>	1	0	1	0
<i>Proteus mirabilis</i>	2	0	2	0
<i>Proteus vulgaris</i>	3	0	3	0
Salmonella spp	11	0	11	0
Salmonella (Group B)	1	0	1	0
Salmonella (Group C)	1	0	1	0
<i>Salmonella flexneri</i>	1	0	1	0
<i>Shigella dysenteriae</i>	1	0	1	0
<i>Shigella sonnei</i>	1	0	1	0
<i>Aeromonas hydrophila</i>	1	0	1	0
<i>Staphylococcus aureus</i>	1	0	1	0
<i>Enterobacter agglomerans</i>	2	0	2	0
<i>Yersinia enterocolitica</i>	1	0	1	0
<i>Clostridium difficile</i>	1	0	1	0
Vibrio spp	3	0	3	0
<i>Vibrio cholera</i>	1	0	1	0
<i>Vibrio parahaemolyticus</i>	1	0	1	0
<i>Hafnia alvei</i>	1	0	1	0
Acinetobacter spp	1	0	1	0
<i>Acinetobacter baumannii</i>	1	0	0	1 (strong)
Yeast	1	0	1	0

All thermophilic Campylobacters tested were reactive with Microgen® Campylobacter. Non-thermophilic campylobacters were non-reactive although 3 of the 11 strains tested showed weak non-specific agglutination.

Of a group of closely related bacteria, 2 *Helicobacter pylori* isolates (of 12) were weakly reactive in the test and 2 isolates of *Helicobacter cinaedi* cross-reacted.

A wide range of other non-related bacteria was also tested. None were reactive in the test although one isolate of *Acinetobacter baumannii* exhibited strong non-specific agglutination.

## REPRODUCIBILITY

**Intra-batch reproducibility** was established by testing sensitivity of 1 batch of product against serial dilutions of reference and kit positive control antigens on 3 separate occasions. End-point titres obtained with reference/control antigens were identical in the three assays.

**Inter-batch reproducibility** was examined by testing sensitivity of 3 batches of product against serial dilutions of reference and kit positive control antigens. No significant differences in end-point titres were seen across the three batches.

## APPENDIX

### Preparation of Blood-Free Selective Medium - Modified CCDA Preston

1. Dissolve 22.75g Campylobacter Blood-Free Selective Agar Base (Oxoid CM739) in 500mL distilled water by boiling.
2. Sterilise by autoclaving at 121°C for 15 minutes.
3. Cool to 50°C, then add 1 vial Cefoperazone Selective Supplement (Oxoid SR125) aseptically.
4. Mix thoroughly and pour into sterile petri dishes (15-20mL/dish)

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Microgen Bioproducts Ltd