

KAP-RPLA "SEIKEN"

A Reagent Kit for detecting *Vibrio parahaemolyticus* heat-stable hemolytic toxin by Reversed Passive Latex Agglutination 20 tests

Demonstration of the producibility of heat-stable hemolytic toxin is regarded as a decisive proof of the human pathogenicity of *Vibrio parahaemolyticus*, and the hemolytic ring occurring during the bacterial growth on Wagatsuma medium which suggests the presence of the toxin is known as the "Kanagawa phenomenon."

The sensitized latex in this test kit is prepared as follows. Rabbits are immunized with the specifically purified heat-stable hemolytic toxin, the antisera taken from the rabbits are purified by affinity-chromatography and polystyrene latex particles are sensitized with the purified antisera.

CHARACTERISTICS

1. This test does not require animals or cell cultures and shows good sensitivity and specificity.
2. Detection sensitivity is superior to the "Kanagawa phenomenon" on Wagatsuma medium.
3. Since it uses the microtiter technique this procedure is easy.

KIT COMPONENTS

1. **Sensitized latex:** 5 ml 1 vial
A suspension of latex particles sensitized with specific antibodies (rabbit IgG) to heat-stable hemolytic toxin. The reagent contains sodium azide as a preservative at a concentration of 0.1 w/v%.
2. **Control latex:** 5 ml 1 vial
A suspension of latex particles sensitized with normal rabbit IgG. The reagent contains sodium azide as a preservative at a concentration of 0.1 w/v%.
3. **Control heat-stable hemolytic toxin (lyophilized):** 0.5 ml equivalent 1 vial
This reagent is lyophilized heat-stable hemolytic toxin, and the agglutination titer is shown on the vial label.
4. **Diluent:** 50 ml 1 vial
Phosphate-buffered saline (PBS) containing 0.5 w/v% bovine serum albumin and the reagent contains sodium azide as a preservative at a concentration of 0.1 w/v%.

INTENDED USE

This reagent kit is intended for use in detecting heat-stable hemolytic toxin produced by *Vibrio parahaemolyticus* grown in the intestine.

PRINCIPLE

Latex particles sensitized with antibodies to *Vibrio parahaemolyticus* toxin react specifically with the toxin produced by the organism in the sample and form an agglutination. This reagent kit employs Reversed Passive Latex Agglutination, based

upon this principle.

PROCEDURES

1. Materials necessary for the test

- (1) **Sample** Culture fluid of *V. parahaemolyticus*
- (2) **Equipment** Microtiter plate (V-type)
Dropper (25 μ l)
Diluter (25 μ l)
Mixer for microtiter plate
Others (black paper, plate cover, marking tape, moisture box)
- (3) **Reagent** KAP-RPLA "SEIKEN"
Mannit Peptone solution containing sodium chloride at a concentration of 5%.

2. Preparation of reagent and sample

1) Reagent preparation

Use sensitized latex, control latex and diluent as supplied. Control enterotoxin should be used after reconstituting with 0.5 ml diluent.

2) Sample preparation

Culture the test bacteria in Mannit-Peptone solution (see below) with sodium chloride added by 5% at 37°C for 18 hours. Centrifuge the culture fluid at 3000 rpm for 20 minutes, and use the supernatant as the test specimen.

Preparation of Mannit-peptone solution containing 5% sodium salt.

Polypeptone	2 g
D-Mannit	0.5 g
NaCl	5 g
Distilled water	100 ml


pH 7.8

3. Reversed Passive Latex Agglutination

- 1) Use two rows of microtiter wells for each specimen. Using a dropper add 25 μ l of diluent to all wells except the first well in each row.
- 2) Drop 25 μ l of specimen into the first well in each row.
- 3) Using two diluters, take 25 μ l of sample and carry out two-fold serial dilutions from the second wells up to the second to last wells.
- 4) Add 25 μ l sensitized latex suspension to each well in one row, and add 25 μ l control latex to all the wells in the other row.
- 5) As a positive control, prepare one well as follows. Add 25 μ l reconstituted control heat-stable hemolytic toxin and sensitized latex to one well in the microtiter plate.
- 6) Shake the microtiter plate well with a microtiter plate mixer
- 7) To avoid evaporation of the solutions in the wells, cover the microtiter plate or place it in a moisture box, then leave it for 18–20 hours and observe the results.
- 8) Observation of the microtiter plate should be carried out over black paper placed in a well-lit place looking from above.

INTERPRETATION

Refer to the following figures in the observation of results.



Regard agglutination stronger than (+) as positive.

If the control latex produces agglutination stronger than (+), then this should be regarded as non-specific agglutination.

If agglutination stronger than (+) is observed in the last well where only latex reagents and diluent were added, agglutination has been spontaneous and the kit should not be used.

It is recommended that the test strain which is interpreted as false positive at lower dilutions be tested after incubation with shaking at 37 °C for 18 hours.

PERFORMANCE CHARACTERISTICS

1. Sensitivity and detection range

Specimen	Dilution							
	× 2	× 4	× 8	× 16	× 32	× 64	× 128	× 256
KAP-T (100 ng/ml)	3+	3+	3+	3+	3+	2+	+	—
	3+	3+	3+	3+	3+	2+	±	—
	3+	3+	3+	3+	3+	2+	+	—
KAP-T 2-fold dilution	3+	3+	3+	2+	2+	+	—	
	3+	3+	3+	2+	+	±	—	
	3+	3+	3+	2+	2+	±	—	
KAP-T 10-fold dilution	3+	2+	+	—				
	3+	2+	+	—				
	3+	2+	+	—				

Marks represent the degree of agglutination.

(in-house data)

When samples of purified *V. parahaemolyticus* toxin (KAP-T, 100 ng/ml) prepared by Denka Seiken were tested using this kit according to the instructions, the undiluted toxin solution, a 2-fold dilution, and a 10-fold dilution produced agglutination up to a 64–128× dilution, 32–64× dilution, and 8× dilution, respectively. All control latex produced negative results, and it was concluded that the sensitized latex had a sensitivity of 1–2 ng/ml.

2. Specificity and within-run reproducibility

Sample	Test									
	1	2	3	4	5	6	7	8	9	10
Positive (1)	× 64	× 64	× 128	× 128	× 128	× 128	× 128	× 128	× 128	× 128
Positive (2)	× 1024	× 512	× 1024	× 1024	× 1024	× 1024	× 512	× 1024	× 1024	× 512
Positive (3)	× 128	× 128	× 128	× 64	× 128	× 64	× 128	× 128	× 128	× 128
Negative (1)	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2
Negative (2)	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2
Negative (3)	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2	< × 2

(in-house data)

Figures in the table represent the maximum specimen dilution factor which produced agglutination.

Culture media of the *Vibrio parahaemolyticus* strains which show positive signs of Kanagawa phenomenon (positive samples 1, 2 and 3) and a culture media of the strain which yields a negative result of this phenomenon (negative samples 1, 2 and 3) were tested with this kit 10 times. All positive samples produced positive results, and all negative samples produced negative results

The maximum dilution factor of the positive reaction did not vary by more than one well in 10.

Control latex produced negative results.

NOTE

1. Use this test kit only for in vitro diagnostic purposes.
2. Do not freeze the kit. Ensure the kit has reached room temperature at least 30 minutes before use.
3. Use a microtiter plate which is free from scratches and stains.
4. Thoroughly agitate the latex reagent vial to resuspend the latex particles and form a homogeneous solution.
5. Store the reconstituted control enterotoxin at 2–10°C . In this manner, the reconstituted reagent can be used for up to 3 months.
6. Do not mix reagents of different lot numbers.
7. Materials and equipment used in this test should be sterilized by soaking in 0.1% sodium hypochlorite solution (chloride content about 1000 ppm) for more than 1 hour or by autoclaving them at 121 °C for more than 20 minutes.
8. Sodium azide contained in the reagents may react with lead and copper piping to form highly explosive metal azides. Upon disposal, flush this reagent away with a large volume of water to prevent accumulation of azides.

STORAGE AND EXPIRATION

Store the reagent at 2–10°C, up to expiration date on the label.

PACKAGE

KAP-RPLA “SEIKEN” 20 tests per box

REFERENCES

1. Miyamoto, Y., et al.: In Vitro Hemolytic Characteristic of *Vibrio parahaemolyticus*: Its Close Correlation with Human Pathogenicity, *J. Bacteriol.*, **100**, 1147 (1969).
2. Sakurai, J., et al.: Purification and Characterization of Thermostable Direct Hemolysin of *Vibrio parahaemolyticus*, *Infect. Immun.*, **8**, 775 (1973).